



MagnaWand™

Cleaning Instruments

A non-magnetic tool designed specifically to clean and disinfect MRI and CAT scan machines. Finally Hospitals will be able to easily clean between patients without losing valuable scanning time.

The **MagnaWand** tool reaches inside the MRI magnet bore or CAT scan gantry to clean and disinfect those difficult to reach places in seconds.


Magnawand's specifically designed disposable pads can be used with any spray disinfectant or cleaning solution.

Once the area of interest has been cleaned and disinfected, simply eject the disposable pad without touching it with your hands, eliminating the risk of further contamination.



Patent Pending

The MagnaWand weighs less than 1 lb and is approx. 36 inches long.



The MagnaWand Acrylic disposable dispenser is self adhesive, can be mounted or stand alone atop a counter. Easily dispenses one disposable at a time. Comes with 30 disposables.

Our Wand hook, is designed to hold our MagnaWand, it is self adhesive and can be mounted on the MRI magnet, or on the side of a counter, next to our Acrylic dispenser. Together they will save you time and space!!!!

Did You Know...

- More than 50 Million MRI and CAT scan studies performed every year
- The Spread of resistant bacteria in Hospitals HAS TRIPLED SINCE 1980
- MRI and CAT scan machines are NOT routinely cleaned between patients
- Deadly bacteria can survive on hard surfaces for up to an hour



Each MagnaWand cleaning pad measures 4" X 5" inches.